

SEREN

Powerchair

*“Available with a
Captain or Clinical seat”*



CE



Crash Tested

The nimble Seren is an ultra-maneuvrable powerchair which can be used both indoors and outdoors. Compact and agile, the Seren moves effortlessly around tight corners, so it's ideal in areas where space is limited.

drive[™]
www.drivemedical.co.uk

The Seren Powerchair uses a dynamic control system to give a superb non-stop driving performance, and features both front and rear suspension to ensure a smooth and comfortable ride. It is available with either a comfortable Captains seat with two tone upholstery, or a fully adjustable seat with tension adjustable backrest and half folding back.



Front and rear suspension to ensure a smooth and safe ride



Dynamic Shark Controller



Clinical / Adjustable Seat



Captain's Seat

Superb Comfort, Style & Performance...

- » Maximum speed of 4mph*.
- » Maximum range of 12 miles on a full battery charge*.
- » Available in 41cm (16")** and 46cm (18") seat widths.
- » Padded Captain's seat and adjustable headrest provide additional comfort and support.
- » Fold down back mechanism for easy transportation and storage.
- » Designed for indoor use and limited outdoor use such as shopping centres.
- » Compact, agile and ultra manoeuvrable for use in tight spaces and busy environments.
- » Front-wheel drive system.
- » Height and width adjustable armrests with PU arm pads to suit comfort of user.
- » Flip up armrests facilitate side transfers.
- » Tension adjustable backrest to suit comfort preferences**.
- » Width, height and depth adjustable seat with cushion and adjustable back width to suit various user sizes**.
- » Front and rear suspension to ensure a smooth and safe ride.
- » Swivel seat facilitates transfers on and off the chair.
- » Aluminium 25cm (10") middle wheels and 15cm (6") front and back castors with solid puncture proof tyres.
- » Single centre mount footplate.
- » Dynamic Shark controller.
- » Powerful 420W motors.
- » Off board charging (2 x 36Ah batteries).
- » Fitted with a seat belt for extra security.
- » This powerchair has been crash tested providing added safety when being transported in a vehicle.
- » 1 year parts warranty and 2 year warranty on frame.
- » Maximum weight capacity 135kg (21st).

* Subject to conditions of use.
** Clinical / Adjustable Seat only.

Additional options available -

- » Attendant Control.
- » Dual Control.
- » Swing Away Joy Stick Holder.
- » Elevating Legrest.
- » Standard Legrest.

SEREN

CODE: SERENBE

Max Weight Capacity:
135kg (21st)



Length:
100cm (39")



Height:
105cm (41")



Width:
61cm (24")



Seat Width:
41/46cm (16/18")



Seat Depth:
45/47/50cm
(18/19/20")



Weight:
80kg (176lb)

Leaflet Code: LL205 If you would like to view this product please contact your local stockist, details below.