
# Drive Medical Aluminium Wheelchair Range











# Our Wheelchair Range at a Glance

#### Travelite Transport Chair



**Expedition Travel Chair** 



Lightweight Aluminium



XS



SD<sub>2</sub>



K-Chair



Travel Chair in a Bag



**Expedition Lite Transit Chair** 



Ultra Lightweight Aluminium



XS<sub>2</sub>



Spirit



## Travelite Transport Chair





## "No other wheelchair folds auite like the TraveLite!"



# Sturdy yet lightweight, this occasional use chair is ideal for many situations from shopping trips to holidays.

- » Attractive and lightweight aluminium frame in gunmetal grey
- » Folds to approx 1/3 the size of a typical folded wheelchair for easy transportation and storage
- Ideal for holidays, days out and shopping trips
- » Flip up padded desk arms to facilitate easy side transfers and provide additional comfort
- Multiple cross braces for added strength and stability
- » 20cm (8") solid puncture proof tyres provide durability and low maintenance
- » Lightweight travel bag with wheels to make transportation easy
- Excellent retail packaging

## Travel Chair In a Bag



## **Expedition Travel Chair**





## "Comes complete with bag as standard"

# The Enigma Travel Chair is lightweight at only 9kg (19lb) making it ideal for transportation.

- » Available in a Metallic Silver finish
- » Padded back canvas with phone and map pockets for storage
- » A half folding back allows easy transport and storage
- Padded desk style armrests provide additional comfort and allow unhindered access to work surfaces
- » Nylon seat canvas is durable and easy to clean
- » Lap belt for added safety
- » Parking brakes for added safety

- » Fixed 20cm (8") puncture proof rear wheels are durable and low maintenance
- » 15cm (6") front castors with puncture proof tyres
- » Quick release swing out footrests assist with transfers on and off the chair
- » Flip up stepper tube allows attendant to safely mount kerbs
- » Compact design with cranked push handles to reduce storage space
- » Lightweight weighing only 19lb (9kg) including all components (exc. bag)

## "Vivid Blue paint finish"

The Expedition chair is the latest refinement in the Drive range of travel chairs. This model features stylish side panels and tool free height adjustable footplates.

- Attractive Silver finish with a lightweight aluminium frame for easy manoeuvrability
- » Half folding backrest allows easy storage and transportation
- » Full length armrests with upholstered arm pads increase user comfort and are easy to wipe clean
- » Nylon 48cm (19") seat and back offer comfortable support and are durable and easy to clean
- 20cm (8") puncture proof tyres are low maintenance
- » Quick release swing away footrests to facilitate transfers on and off the chair
- Lap belt for added safety
- » Tool free footplate height adjustment
- » Two stepper tubes enable attendant to safely raise front wheels for mounting kerbs

## **Expedition Lite Transit Chair**



## Lightweight Aluminium





### The Expedition Plus is our most feature rich transit chair providing excellent practicality and comfort.

- Attractive Silver finish with a lightweight aluminium frame for easy manoeuvrability
- Arthritic friendly cable brakes allow immediate access for attendant to reduce speed or lock brakes for parking
- Additional parking brakes for extra safety
- Half folding backrest allows easy storage and transportation
- Full length armrests with upholstered arm pads increase user comfort and are easy to wipe clean
- Nylon seat and higher backrest offer comfortable support and are durable and easy to clean
- Larger 30cm (12") puncture proof rear solid tyres are durable and low maintenance
- Quick release swing away footrests to facilitate transfers on and off the chair
- Lap belt for added safety
- Two stepper tubes enable attendant to safely raise front wheels for mounting kerbs
- Lightweight, weighing only 11.9kg (including footrests)

The Lightweight Aluminium Wheelchair offers a value for money introduction to the Aluminium Wheelchair range without compromising on quality. Lightweight and foldable with a half folding back, quick release rear wheels and footrests.

- Constructed from sturdy aluminium, this wheelchair is durable and robust, yet lightweight for transporting (weighing only 9kg (19lb) in its carrying format)
- Nylon seat and back offer comfortable support and are easy to clean
- Chair folds up easily for storage and transportation, further enhanced by half folding backrest, quick release footrests and rear wheels (self propelling models only)
- Upholstered arm pads increase user comfort and are easy to wipe clean
- Tool free height adjutable footrest to 5 positions to optimise user comfort and support
- Footrests can also be swung out of the way to facilitate transfers

- Rear stepper tubes allow carers/ attendants to safely climb kerbs
- Puncture proof front wheels ensure low maintenance and longevity
- Supplied with lap belt for increased user safety

#### Self Propel Model

- » Aluminium 61cm (24") rear wheels with handrims and grey pneumatic tyres allow a smooth independent ride for active users whilst making the chair lightweight
- Total weight: only 15kg (32lb)

- Cable brakes on rear handles allow immediate access for attendant
- Total weight: only 12kg (27lb) for easy transportation

## Ultra Lightweight Aluminium



# "This ultra lightweight wheelchair features a stylish new design"

# The Ultra Lightweight Wheelchair combines style with function and practicality.

- Strong ultra lightweight frame is stable, yet very lightweight for transporting and easily folded for storage
- Half folding backrest, quick release footrests and rear wheels (self propelling model only) further facilitates storage and transportation
- » Grey two tone upholstery offers comfortable support, is easy to clean and along with Silver frame has an attractive appearance
- » Soft fabric side panels improve comfort
- » Features handy rear pocket for storing items
- » Footrests easily adjust into one of five positions to suit user support requirements and swing out of the way to facilitate transfers
- » Desk style armrests allow easy access to work surfaces
- » Two stepper tubes enable attendant to safely raise front wheels for kerbs

- » Two solid Polyurethane (PU) front tyres, 20cm (8") in diameter, are durable, low maintenance and assist with manoeuvrability
- » Long reach brakes on rear tyres are easily activated by user or attendant for parking

#### Self Propel Model

» Spoked 61cm (24") rear wheels with hard wearing grey Polyurethane (PU) tyres and low maintenance, coupled with light aluminium rims and hand rims increase manoeuvrability

#### Transit Mode

- » Cranked push handles allow compact folding
- » Cable brakes on rear handles provide attendant with immediate access
- » Durable 32cm (12.5") solid polyurethane (PU) rear wheels assist with manoeuvrability

## **XS Aluminium**



## "Available in the stylish Spirit Red colour"



# The XS Aluminium Wheelchair has many transportation friendly, quick release components and built in comfort features.

- » Choice of 61cm(24") quick release wheels (self propel) or 32cm(12.5") fixed wheel (transit)
- » Sturdy, colourful aluminium frame is both lightweight for transporting and durable
- » Simple folding mechanism, half folding backrest and removable arms as well as quick release footrests and self propelling rear wheels allow space saving storage
- » Footrests swing in and out to enable transfers on and off chair, whilst flip up armrests facilitate side transfers
- » Desk style armrests provide additional comfort and allow easy access to work surfaces

- Seat and back cushions with tension adjustment for comfort
- Rear stepper tubes with integrated height adjustable anti tip wheels facilitate safe kerb climbing for both attendant and independent user
- » Rear wheels can be adjusted into different position
- » Solid puncture proof Polyurethane (PU) tyres offer durability and minimum maintenance
- » Available in two different seat widths to suit various user sizes
- Lightweight frame weighing under 9kg (20lb) in carry format

#### Colour Options



Spirit Red



Blue

## XS<sub>2</sub> Aluminium





## "Available in 3 seat widths"

# Our new and improved XS2 Wheelchair uses a similar frame to the XS2 Wheelchair but adds many new features.

- » Sturdy, Silver aluminium frame is lightweight for transporting and durable
- Simple folding mechanism, half folding backrest, removable arms as well as quick release footrests and self propelling rear wheels allow space saving storage and easy transportation
- » Adjustable desk style armrests provide additional comfort, flexibility and allow unhindered access to work surfaces
- » Footrests swing in and out to enable transfers on and off chair, whilst flip up armrests facilitate side transfers
- Incorporated as standard are a seat cushion, tension adjustable backrest and push button length adjustment on armrests, allowing chair to be easily adapted to suit comfort preferences

- » A one touch lap belt is also supplied on both self propel and transit models, providing additional stability for user
- » Long reach brakes on rear wheels give easy access to attendant and user for increased safety, with transit model also featuring 'immediate access' cable brakes to reduce speed whilst moving
- » Rear stepper tubes with integrated anti tip wheels facilitate safe kerb climbing for both attendant and independent user
- » Chair is fitted with reflectors to improve outdoor visibility at dusk or when dark
- » All round solid Polyurethane (PU) tyres offer durability and minimum maintenance
- » Available in three different seat widths to suit various user sizes

## **SD2** Aluminium



## "For those who require maximum features"

# The SD2 Wheelchair is an upgrade in both design and features and replaces the old Drive Super Deluxe Wheelchair.

- » Elegant metallic black aluminium frame
- » Double cross brace provides added support to the wheelchair
- » Attendant or user accessible parking brakes on rear wheels ensure safety along with attendant cable brakes on transit models
- Improved footrests swing in and out to enable transfers on and off chair, whilst flip up and removable armrests facilitate side transfers
- Push button adjustment facilitates height and length adjustable arm pads
- » Half folding back mechanism allows easy storage and transportation
- » Incorporated as standard are seat and back cushions and a padded tension adjustable backrest for added comfort

- » Rear stepper tubes with integrated height adjustable anti tip wheels facilitate safe kerb climbing for both attendant and independent user
- » Quick release rear wheels on Self Propel models
- » Black pneumatic rear tyres and puncture proof front tyres are durable and low maintenance
- » Chair is fitted with reflectors to improve outdoor visibility at dusk or when dark
- One touch lap belt for added safety
- » Features handy rear pocket for storing items

## **Spirit**



## K-Chair



## "The ultimate in style and design"

The Spirit Wheelchair uses the latest in wheelchair design to produce an elegant yet highly practical wheelchair. The sculptured frame is constructed from lightweight aluminium for durability and responsiveness.

- » Attractive and durable frame available in various colours
- » One touch adjustable armrests for easy transfers and comfort
- » Folding frame with removable footrests and half folding back for easy storage
- » Height adjustable footrests swing in and out to enable transfers on and off the chair, whilst one touch adjustable armrests facilitate side transfers
- » Stylish 61cm (24") mag wheels with strong, resilient and lightweight plastic handrims
- » Incorporated as standard are padded nylon seat and back cushions

- and a tension adjustable backrest for added comfort
- » Solid puncture proof front and rear tyres provide durability and low maintenance
- » Attendant or user accessible parking brakes on rear wheels ensure safety
- » Rear stepper tubes with integrated height adjustable anti tip wheels facilitate safe kerb climbing for both attendant and independent user
- » Quick release rear wheels for easy storage and transportation
- » Available with mag wheels (Spirit Red only)

## "Adjustable suspension for optimum comfort"

The K-Chair has been developed to provide a new level of comfort for wheelchair users. It features adjustable suspension that can be tailored for the user and the environment thus reducing the effects of any uneven floor surface.

- » Folding frame with removable footrests and half folding back mechanism for easy storage and transportation
- » Padded desk style armrests provide additional comfort and allow unhindered access to work surfaces
- » Incorporated as standard are padded seat and back cushions and a tension adjustable backrest for added comfort
- » Shock absorbing front castors and rear adjustable shock absorbers for a smooth ride and added comfort
- » 61cm (24") rear wheels (quick release) with hand rims, black puncture proof solid tyres and 6 position axle

- Height adjustable footrests with new easy to use pull up handle swing easily in and out to enable transfers on and off the chair
- » One touch adjustable armrests facilitate side transfers
- » Height adjustable anti tip wheels facilitate safe kerb climbing for both the attendant and independent user
- » Stylish ergonomic hand grips enable safe handling for the attendant
- » Attendant or user accessible parking brakes on rear wheels ensure safety

#### Colour Options



Spirit Red



Electric Blue



Silver





Spirit Red



Midnight Blue

# SPECIFICATION

CODE	DESCRIPTION	WIDTH (OPEN)	WIDTH (CLOSED)	LENGTH (INC FOOTRESTS)	HEIGHT	SEAT WIDTH	SEAT LENGTH	SEAT HEIGHT	WEIGHT (INC FOOTRESTS)	WEIGHT CAPACITY
TCoo5	Travelite Transport Chair 18"	58cm (23")	26cm (10")	97cm (38")	93cm (37")	46cm (18")	42cm (16.5")	48cm (19")	9.2kg (20lb)	115kg (18st)
TCoo2SIL	Lightweight Travel Chair In a Bag 19"	58cm (23")	22cm (9")	104cm (41")	90cm (35")	48cm (19")	39cm (15")	46cm (18")	8.5kg (191b)	115kg (18st)
EXPoo <sub>1</sub> SIL	Expedition Travel Chair 18"	54cm (21")	27cm (11")	95cm (37")	97cm (39")	46cm (18")	40cm (16")	46cm (18")	10kg (22lb)	115kg (18st)
EXPoo2SIL	Expedition Lite Transit Chair 18"	54cm (21")	23cm (9")	95cm (37")	96cm (38")	46cm (18")	40cm (16")	50cm (20")	11.9kg (26lb)	135kg (21st)
LAWCoo1	Lightweight Aluminium Wheelchair 18" Self Propel	64cm (25")	34cm (13")	106cm (42")	90cm (35")	46cm (18")	40cm (16")	48cm (19")	14.4kg (32lb)	115kg (18st)
LAWC002	Lightweight Aluminium Wheelchair 18" Transit	56cm (22")	28cm (11")	96cm (38")	92cm (36")	46cm (18")	40cm (16")	48cm (19")	12.3kg (27lb)	115kg (18st)
LAWC007A	Ultra Lightweight Aluminium 18" Self Propel	63cm (25")	29cm (11")	107cm (42")	93cm (37")	42cm (17")	40cm (16")	50cm (20")	14.7kg (32lb)	115kg (18st)
LAWC008A	Ultra Lightweight Aluminium 18" Transit	61cm (24")	28cm (11")	104cm (41")	96cm (38")	42cm (17")	39cm (16")	49cm (20")	12.8kg (28lb)	115kg (18st)
LAWConA	Ultra Lightweight Aluminium 20" Self-Propel	70cm (28")	29cm (11")	107cm (42")	93cm (37")	50cm (20")	40cm (16")	50cm (20")	15.3kg (34lb)	125kg (20st)
LAWC012A	Ultra Lightweight Aluminium 20" Transit	66cm (26")	28cm (11")	104cm (41")	96cm (38")	50cm (20")	39cm (16")	49cm (20")	13.4kg (29lb)	125kg (20st)
XSAWCTS18ST	XS Aluminium 18" Transit	63cm (25")	26cm (10")	104cm (41")	92cm (36")	46cm (18")	43cm (17")	52cm (20")	15kg (32lb)	115kg (18st)
XSAWCSP18ST	XS Aluminium 18" Self Propel	62cm (25")	33cm (13")	104cm (41")	94cm (37")	45cm (18")	40cm (16")	52cm (20")	16.7kg (37lb)	115kg (18st)
XSAWCTS20ST	XS Aluminium 20" Transit	68cm (27")	26cm (10")	104cm (41")	92cm (36")	50cm (20")	43cm (17")	52cm (20")	16.7kg (37lb)	115kg (18st)
XSAWCSP2oST	XS Aluminium 20" Self Propel	67cm (27")	33cm (13")	104cm (41")	94cm (37")	50cm (20")	40cm (16")	52cm (20")	16.9kg (37lb)	115kg (18st)
XS2SP16SIL	XS2 Aluminium 16" Self-Propel	58cm (23")	30cm (12")	103cm (40.5")	93cm (37")	40cm (16")	43cm (17")	52cm (20.5")	16.4kg (36lb)	135kg (21st)
XS2TS16SIL	XS2 Aluminium 16" Transit	58cm (23")	29cm (11")	102cm (40")	95cm (37")	40cm (16")	43cm (17")	52cm (20.5")	15.6kg (34lb)	135kg (21st)

XS2SP18SIL	XS2 Aluminium 18" Self-Propel	63cm (25")	30cm (12")	102cm (40")	93cm (37")	45cm (18")	43cm (17")	52cm (20.5")	16.6kg (37lb)	135kg (21st)
XS2TS18SIL	XS2 Aluminium 18" Transit	63cm (25")	29cm (11")	102cm (40")	95cm (37")	45cm (18")	43cm (17")	52cm (20.5")	15.8kg (34lb)	135kg (21st)
XS2SP2oSIL	XS2 Aluminium 20" Self-Propel	68cm (27")	30cm (12")	102cm (40")	93cm (37")	50cm (20")	43cm (17")	52cm (20.5")	16.8kg (37lb)	135kg (21st)
XS2TS20SIL	XS2 Aluminium 20" Transit	68cm (27")	29cm (11")	102cm (40")	95cm (37")	50cm (20")	43cm (17")	52cm (20.5")	16kg (34lb)	135kg (21st)
SD <sub>2</sub> SP <sub>1</sub> 6BLK	Super Deluxe 2 16" Self Propel	61cm (24")	32cm (12.5")	105cm (41")	92cm (36")	40cm (16")	40cm (16")	49cm (20")	19.8kg (44lb)	135kg (21st)
SD2TS16BLK	Super Deluxe 2 16" Transit	57cm (22")	28cm (11")	97cm (38")	96cm (38")	40cm (16")	40cm (16")	49cm (20")	16.4kg (36lb)	135kg (21st)
SD <sub>2</sub> SP <sub>1</sub> 8BLK	Super Deluxe 2 18" Self Propel	66cm (26")	32cm (12.5")	105cm (41")	92cm (36")	45cm (18")	40cm (16")	49cm (20")	20kg (44lb)	135kg (21st)
SD <sub>2</sub> TS <sub>1</sub> 8BLK	Super Deluxe 2 18" Transit	62cm (25")	28cm (11")	97cm (38")	96cm (38")	45cm (18")	40cm (16")	49cm (20")	16.6kg (37lb)	135kg (21st)
SD2SP20BLK	Super Deluxe 2 20" Self Propel	71cm (28")	32cm (12.5")	105cm (41")	92cm (36")	50cm (20")	40cm (16")	49cm (20")	20.2kg (44lb)	135kg (21st)
SD2TS20BLK	Super Deluxe 2 20" Transi	67cm (26")	28cm (11")	97cm (38")	96cm (38")	50cm (20")	40cm (16")	49cm (20")	16.8kg (37lb)	135kg (21st)
SD2SP22BLK	Super Deluxe 2 22" Self Propel	76cm (30")	32cm (12.5")	105cm (41")	92cm (36")	55cm (22")	40cm (16")	49cm (20")	20.4kg (44lb)	135kg (21st)
SD2TS22BLK	Super Deluxe 2 22" Transit	72cm (28")	28cm (11")	97cm (38")	96cm (38")	55cm (22")	40cm (16")	49cm (20")	17kg (37lb)	135kg (21st)
XSES18RD	Spirit Wheelchair 18" - Spirit Red	71cm (28")	32cm (13")	102cm (40")	91cm (36")	45cm (18")	40cm (16")	48cm (19")	17.3kg (38lb)	115kg (18st)
XSES18BL	Spirit Wheelchair 18" - Electric Blue	71cm (28")	32cm (13")	102cm (40")	91cm (36")	45cm (18")	40cm (16")	48cm (19")	17.3kg (38lb)	115kg (18st)
XSES18SI	Spirit Wheelchair 18" - Silver	71cm (28")	32cm (13")	102cm (40")	91cm (36")	45cm (18")	40cm (16")	48cm (19")	17.3kg (38lb)	115kg (18st)
XSESMAG18RD	Spirit Wheelchair 18" with MAG Wheels - Red	71cm (28")	32cm (13")	102cm (40")	91cm (36")	45cm (18")	40cm (16")	48cm (19")	17.3kg (38lb)	115kg (18st)
XSSW18RD	Enigma K-Chair 18" - Spirit Red	66cm (26")	35cm (14")	107cm (42")	98cm (39")	45cm (18")	42.5cm (17")	48cm (19")	19kg (42lb)	115kg (18st)
XSSW18BL	Enigma K-Chair 18" - Midnight Blue	66cm (26")	35cm (14")	107cm (42")	98cm (39")	45cm (18")	42.5cm (17")	48cm (19")	19kg (42lb)	115kg (18st)

## Comfort, Style, Performance & Reliability.

We understand the true value of innovation within this market place; and know that our Wheelchair development will have a direct influence on improving an individual's life.

Leaflet Code: LL137

(€

All sizes and weights stated are nominal. All errors and omissions excepted.

Due to the limitations of photography and the printing process colours shown may not be 100% accurate.

