

Lightweight paediatric wheelchair

Equipped with height adjustable push handles Easy foldable and very compact





Very easy to fold, very compact and easily transported.



Easy swing inwards/outwards removable footrests.



9 kg

Including skirtguard and optional skirtguard with spray guard.



Height adjustable push handles.

Standard specifications

Seat width: 30 and 35 cm Seat depth: 30 and 35 cm 401/2 and 43 cm Seat height: Backrest height: 38, 40 and 42 cm Total width: Seat width + 21 cm Total length: > 83 cm Total height: > 80 cm Width folded: 28 cm

20 cm

6" x 1" (other sizes available)
22" x 13%" (other sizes available) Front castor size: Rear wheel (sp):

Armrest height:

Total weight: 13.9 kg Weight (transport): Max. user weight: 75 kg

> Frame colour available in:



Product complies with crash test ISO 7176-19:2001 and EN12183



Corunna House 42-44 Ousegate Selby Y08 4NH

Van Os Medical UK Ltd

Tel.: +44 (0) 1757 701177 Fax: +44 (0) 1757 706011

sales@vanosmedical.com