

Lightweight paediatric wheelchair

Multiple configurations from a single frame, which can be adjusted as the child grows.



Available in metallic blue, metallic pink and metallic green

Standard specifications

Seat width: 30, 32½ and 35 cm
Seat depth: 30, 32½ and 35 cm
Seat Height: 40, 42, 44½, 45½, 46 and 46½ cm

Backrest height: 40, 42½ and 45 cm (configurable)
Total width: Seat width + 21 cm
Total length: 86 cm
Total height: 84 cm
Width folded: 33 cm
Armrest height: 13-23 cm

Rear wheel (sp): 20" or 22" x 1½" (depends on seat height)

Total weight (complete): 14.4 kg
Weight (transport): 8.4 kg
Max. user weight: 80 kg

Frame colour available in:



Product complies with crash test ISO 7176-19:2001



Van Os Medical UK Ltd
Excel House
Escrick Business Centre
Escrick, York
YO19 6FD

Tel.: +44 (0) 1904 720170

sales@vanosmedical.com