

 **Breezy**<sup>®</sup>  
Premium products. Quality service.

# *Moonlite*



*Moonlite*



CRASH TESTED



## Compact and lightweight....

Whether you're just taking a quick trip to the shops or going out for the day with your loved ones, the Breezy Moonlite is the ideal partner. Lifting it into the car or storing it in the tightest of spaces is effortless thanks to the revolutionary compact design and lightweight frame.

## Added comfort...

Comfort is always a priority when choosing a wheelchair and with the Breezy Moonlite's carefully designed seating position and padded upholstery, from the moment you sit down you'll know it's the chair for you.

## Peace of mind...

The Breezy Moonlite leads the way in terms of safety, quality and durability. Built to last, the Breezy Moonlite has passed crash testing ensuring you have one of the safest products on the market.

## User friendly features...



Attendant brake on attendant propelled model



Swing-in / swing-out legrests



Half-folding backrest



Folds for storage



Lap belt




# Moonlite




41cm, 46cm  
(16", 18.1")




SW + 19cm (7.4") (AP)  
SW + 21cm (8.2") (SP)




41cm (16")



48cm  
(18.8")



48cm  
(18.8")



96cm (37.7")

## Technical data sheet

<b>Upholstery colour:</b>	<b>black / grey</b>
<b>Approved user weight:</b>	<b>115kg (18 stone)</b>
<b>Chair weight:</b>	<b>11.8kg (AP) 13.5kg (SP)</b>
<b>Folded width: (with wheels)</b>	<b>23.5cm (9.2") (AP) 30cm (11.8") (SP)</b>
<b>Transit tested to ISO 7176 Part 19 and passed:</b>	<b>Yes</b>

AP - Attendant propelled, SP - Self propelled

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.com](http://www.SunriseMedical.com) or alternatively is available on request in large text.



Sunrise Medical Limited  
Thorns Road,  
Brierley Hill  
West Midlands  
DY5 2LD

Tel: +44 (0) 845 605 66 88  
Fax: +44 (0) 845 605 66 89

[www.SunriseMedical.com](http://www.SunriseMedical.com)