

Taking you anywhere, anytime Ergo Lite 2 From only 8.7 kg

## Erqo Lite 2

The KM-2512 is an ultra light wheelchair which weighs from just 8.7kg complete and breaks down to only 6.4kg (self propel version).

It is a durable, lightweight wheelchair which comes with Karma's patented S-Ergo seating system for enhanced comfort. The S-Ergo seat helps redistribute

pressure and also helps stop the user sliding down the seat.

The removable swing in/out footrest hangers are convenient to use and allow the KM-2512 to collapse into a compact package.





Swing in/out detachable footrest hangers



KM-2512 F14



The lowest transportation weight is only 6.4 kg (without footrest hangers and rear wheels)

## **FEATURES**

- S-Ergo system provides pressure relief relief and helps stop the user sliding down the seat
- AEGIS Microbe Shield permanently bonded antimicrobial barrier
- Swing in/out detachable footrest hangers
- Foldable backrest and frame for effortless storage and transportation
- Desk-length padded armrests
- Attendant brakes
- Comfortable & durable upholstery
- 20" Quick release spoke rear wheels

Model	width	F/R Wheels		<u>O</u> T			<u>O</u>	<u>Jr</u>	<u>br</u>		<u>F</u>		Max. User Weight
KM-2512 F14	16"	6" - 14"	41	42	62	27	21	93	48	42	91	8.7 kg	100 kg
KM-2512 F14	18"	6" - 14"	46	42	62	27	21	93	48	42	91	8.8 kg	100 kg
KM-2512 Q20	16"	6" - 20"	41	42	65	31	21	93	48	42	91	10.4 kg	100 kg
KM-2512 Q20	18"	6" - 20"	46	42	65	31	21	93	48	42	91	10.5 kg	100 kg

KARMA reserves the right to modify information herein without further notice



KM-2512 Q20

CE 🖅 🚱 ISO 9001  $\pm 1$  cm  $\pm 0.5$ kg

## **KARMA Mobility Limited**

www.karmamobility.co.uk e-mail:info@karmamobility.co.uk Unit 6, Target Park, Shawbank Road, Redditch, B98 8YN U.K. TEL: 0845-630-3436 FAX: 0845-630-3736