

Invoice to:
 Name
 Street
 Town: Post code:
 E-Mail:
 Tel:

Deliver to:
 Name
 Street
 Town: Post code:
 E-Mail:
 Tel: Fax:

Order date
 Purchase order:

Marked for:

Fax number: +44 (0)845 605 66 89

Email: orders@sunmed.co.uk

Maximum user weight limit 125kg

Base Price	
£835	
Quantity	

STANDARD FEATURES
Standard padded or depth adjustable seat
Tension adjustable backrest upholstery
Height adjustable backrest: 41cm - 46cm (43cm Std)
Flip back removable depth & height adjustable arms
Seat height adjustable from 42.5cm to 50.5cm
6", 7" or 8" castors
12", 22" or 24" rear wheels
Angle and reversible footplates
Seat belt

Crash tested and passed without tie down brackets

CLIENT MEASUREMENT

		Seat Width					
		38cm	40.5cm	43cm	45.5cm	48cm	52cm
Seat Depth	41cm	○*	○*	○*	○	○	○
	43.5cm	○*	○*	○*	○*	○*	○
	46cm	○	○*	○*	○*	○*	○*
	48.5cm [†]	○	○	○	○*	○*	○*
	51cm [†]	○	○	○	○*	○*	○*

* Indicates seat size is available with JAY cushion. Choose below

[†] Seat depth 48.5cm and 51cm only possible with single post height adjustable armrest

Front seat height in cm

Castor & Wheel size	Front seat height in cm												
	37.5	39	40	41	42.5	44	45	46.5	48	49	50.5	52*	53*
6" Castor / 12" rear wheel							○						
7" Castor / 12" rear wheel							○						
8" Castor / 12" rear wheel							○						
6" Castor / 22" rear wheel	○	○	○	○	○	○	○	○	○				
7" Castor / 22" rear wheel				○	○	○	○	○	○	○			
8" Castor / 22" rear wheel					○	○	○	○	○	○	○	○	
5" Castor / 24" rear wheel				○	○	○	○	○					
6" Castor / 24" rear wheel			○	○	○	○	○	○					
7" Castor / 24" rear wheel				○	○	○	○	○	○				
8" Castor / 24" rear wheel					○	○	○	○	○	○	○	○	○

NEW 5" Castor Now Available

* Maximum user weight 125Kg

SEAT UPHOLSTERY			£
RUB020001	○	Standard padded nylon	(Std default option)
RUB020003	□	Nylon depth adjustable	35
Cushions			
RUB020004	□	Cushion 5cm	28
JAY020002	□	JAY Basic	(* available in selected sizes as above) 55
JAY020003	□	JAY Soft Combi P	(* available in selected sizes as above) 77
JAY020001	□	JAY Easy Visco	(* available in selected sizes as above) 132

BACKREST UPHOLSTERY			
RUB030003	<input type="radio"/>	Tension adjustable back (7 strap) <i>(not with J3 backs)</i>	
BACKREST			
RUB030017	<input checked="" type="radio"/>	Fixed push handles	
RUB030019	<input type="radio"/>	Half folding back <i>(not with drum brake)</i>	
RUB030017	<input type="radio"/>	Fixed back - non folding <i>(Std default option)</i>	
RUB030018	<input type="checkbox"/>	Integrated height adjustable push handles <i>(not with half folding back or extension back)</i>	93
RUB030020	<input type="checkbox"/>	Reclining back <i>(Anti-tips recommended; not with seat depth 51cm or rear wheel extension. Only with single post armrest)</i>	523
RUB030021	<input type="checkbox"/>	Angle adjustable back with integrated height adjustable push handles <i>(anti-tips are recommended!)</i>	168
RUB030023	<input type="checkbox"/>	Auto folding stabiliser bar (HeliX style) <i>(not with folding back, headrest; recommended with recliner)</i>	108
HEADREST			
RUB030025	<input type="checkbox"/>	Standard headrest <i>(includes swing-away stabiliser bar)</i>	234
RUB030133	<input type="checkbox"/>	Extension back (domed headrest) <i>(not with int. height adjustable push handles and SW 38-52)</i>	121
ARMREST			
RUB040002	<input type="radio"/>	Flip back std. short & long armpad position (sliding) <i>(not with recliner; not with SD 48.5cm or 51cm)</i>	
RUB040003	<input type="checkbox"/>	Flip Back, height adjustable (T-Post) short & long armpad position (sliding) <i>(not with recliner; not with SD 48.5cm or 51cm)</i>	76
RUB040004	<input type="checkbox"/>	Single post, height adjustable with short & long armpad position (sliding) <i>(compulsory with recliner)</i>	85
RUB040005	<input type="checkbox"/>	Hemi arm pad RelaX 2 style 1 unit (left & right usable) <i>(n/a with tray; only with height adj. armrest)</i>	76
RUB040006	<input type="checkbox"/>	Hemi arm pad RelaX 2 style 2 units (left & right usable) <i>(n/a with tray; only with height adj. armrest)</i>	50
RUB040007	<input type="checkbox"/>	Hemi arm pad large RelaX 2 style 1 unit (left & right usable) <i>(n/a with tray; only with height adj. armrest)</i>	97
RUB040008	<input type="checkbox"/>	Hemi arm pad large RelaX 2 style 2 units (left & right usable) <i>(n/a with tray; only with height adj. armrest)</i>	50
RUB040107	<input type="checkbox"/>	Sideguard, composite, removable, 210 mm high (310 mm long), black	97
LEGREST & FOOTPLATES			
RUB050107	<input checked="" type="radio"/>	70° hanger	
RUB050108	<input type="radio"/>	80° hanger <i>(not in combination with 7" and 8" castors)</i>	
RUB050001	<input type="radio"/>	Angle adjustable footplate	
RUB050022	<input type="checkbox"/>	Footboard, divided, aluminium, angle-adjustable, flip-up (sideways), heel loop, black	97
FMGMD06-1292-02	<input type="checkbox"/>	Divided footplate for short lower leg length	167
RUB050003	<input type="checkbox"/>	Elevating Legrest - PAIR <i>(curved calf pad)</i>	254
RUB050004	<input type="checkbox"/>	Elevating Legrest - LH only	127
RUB050005	<input type="checkbox"/>	Elevating Legrest - RH only	127
RUB050006	<input type="radio"/>	Comfort calf pad (pair)	
RUB050008	<input type="checkbox"/>	Amputee Support - PAIR <i>(Depth, Angle, Height adjustable / Swing-away)</i>	219
RUB050009	<input type="checkbox"/>	Amputee Support - LH only <i>(Depth, Angle, Height adjustable / Swing-away)</i>	121
RUB050010	<input type="checkbox"/>	Amputee Support - RH only <i>(Depth, Angle, Height adjustable / Swing-away)</i>	121
RUB050300	<input type="checkbox"/>	Footboard mounting high mount	140

		£
CASTOR		
RUB080101	○ 8" x 1 1/4" solid	
RUB080102	○ 8" x 50 solid	
RUB080103	○ 8" x 1 1/4 pneumatic	
RUB080104	○ 7" solid	
RUB080105	○ 6" soft solid	
RUB080006	○ 6" solid	
RUB080057	○ 5" soft solid	
REAR WHEELS		
RUB070001	○ 24" x 1" solid (spoked wheels)	
RUB070002	○ 24" x 1" pneumatic (spoked wheels) <i>(not with drum brake)</i>	
RUB070083	□ 24"X 1", pneumatic Right Run, 24 spokes, straight, 125kg max user	270
RUB070003	○ 22" x 1 3/8" solid (spoked wheels) <i>(also for drum brake)</i>	
RUB070004	□ Handrim cover	92
RUB070005	○ 12" solid (fixed axle) <i>(not with drum brake or unilateral wheel lock)</i>	
RUB070006	○ Rear wheel extension 2,5 cm (behind back tube) (3) <i>(n/a with angle adj. or reclining back; 22/24" only)</i>	
RUB070007	□ One arm drive (fixed axle) - LH <i>(pneumatic tyres only)</i>	420
RUB070008	□ One arm drive (fixed axle) - RH <i>(pneumatic tyres only)</i>	420
WHEEL LOCKS		
RUB060001	○ Push to lock	
RUB060002	□ Lever extension (pair)	58
RUB060003	□ Drum brake <i>(24"/22" only spoked wheel solid, not with folding back)</i>	328
RUB060005	□ Unilateral wheel lock - RH	127
RUB060004	□ Unilateral wheel lock - LH	127
ACCESSORIES		
Crutch holder (5)		
RUB090102	□ LH Mounted	48
RUB090103	□ RH Mounted	48
Spoke guards (transparent)		
RUB090004	□ 22" Rear Wheel	122
RUB090105	□ 24" Rear Wheel	122
Swing-away tray (1,2)		
RUB090111	□ LH Mounted <i>(n/a with hemi armrests; all seat widths)</i>	203
RUB090112	□ RH Mounted <i>(n/a with hemi armrests; all seat widths)</i>	203
RUB090110	□ Slide-on tray (2) <i>(n/a with hemi armrests; n/a with 52cm seat width)</i>	121
Anti-tips (4,5)		
RUB090113	□ Pair <i>(pair, recommended with recliner)</i>	52
RUB090114	□ Left <i>(not compatible with step tubes)</i>	26
RUB090115	□ Right <i>(not compatible with step tubes)</i>	26
Step tube		
RUB090117	□ Left <i>(not compatible with anti-tips)</i>	28
RUB090118	□ Right <i>(not compatible with anti-tips)</i>	28
RUB090020	□ Hemi frame - (Front castor outrigger for stability or foot propelling)	70
RUB090109	○ Tool kit	
RUB090019	● Positioning belt (Push button style)	
COLOURS		
RUB090023	○ Brilliant Silver	
RUB090024	○ Matt Blue	
RUB090025	○ Brilliant Red	
RUB090024	○ Matt Blue	
RUB090025	○ Brilliant Red	
NON COMPATIBILITIES		
(1) Only with height adjustable armrests. (2) NA with hemi armrests (3) NA with 48.5 or 51cm seat depth and anti-tips		
(4) Not possible with rear wheel extension and 48.5 or 51cm seat depth (5) Crutch holder with anti-tips not available on seat depth 51cm		
<p>Note: Sunrise Medical strongly recommends the use of a step tube on all transit wheeled models, the XL model and also where attendant push is the predominant intended use.</p>		
Breezy RubiX² - Base Price		£835
Option price		£
Total price		£