

e-goes K5 Folding Powered Chair

The ultimate folding powered chair
for adjustability and features

12 Rear Wheels and Front Suspension



e-goes
SMART E-MOBILITY

SPECIFICATIONS:



Folding Size:
40x62x84 cm / 17.3"x24.4"x33.1"



Unfold Size L*W*H	102x62x96 cm / 40.2"x24.4"x37.8"	Folding Size L*W*H	44x62x84cm / 17.3"x24.4"x33.1"
N.W. (with battery)	35kg/77 lbs	N.W. (without battery)	32.5kg/71.5 lbs
Front Wheel	8" Front Wheel	Rear Wheel	12" Rubber Tyre
Lithium Battery	12Ah	Charging Period	5 Hours
Max Speed	6 KM/H	Motor Power	250W x 2
Cruiser Mileage	15km/9.3mi	Motor Type	Brushless Hub Motor
Brake	Electromagnetic brake	Braking Distance	≤ 1M
Loading Capacity	150 KG	Charger	DC 29.4 V
Turning Radius	960 mm	Climbing Slope	12°

Length	Width	Height	Seat Depth	Seat Width	Seat Height	Backrest Height	N.W.
1020mm 40.2in	620mm 40.2in	960mm 37.8in	430mm 16.9in	460mm 18.1in	520mm 20.5in	500mm 19.7in	32.5kg/71.5lb (without battery)

SMART CONTROLLING SYSTEM

- The LCD Joystick control panel offers several customisable functions including the ability to adjust the responsiveness of the chair from Soft, though Normal to Powerful making the chair easy to set up to suit the conditions and user preferences.

- USB Port



Whole backrest can be reclined
Max 125 degree



Backrest adjustable



Ergonomic armrest to be adjusted up and down, to be adjusted forward and backward .



