

The Dash Stowaway is ideal for users & carers wanting quality, lightness & flexibility not found in traditional wheelchairs



Dash Stowaway *Technical Information*

WEIGHT INFORMATION	
Max occupant weight	110kgs
SEATING AREA DIMENSIONS	
width	17"
depth	17"
backrest canvas height	410mm
backrest frame angle from vertical	10 degrees
seat frame angle from horizontal	5 degrees
height of armrests	213mm
footrest to seat adjustable height	400mm» 470mm
seat front edge to ground	490mm
push handle height	927mm
CORRIDOR WIDTHS TO TURN (315mm wheel)	
through 360 degrees between walls	1425mm
OVERALL DIMENSIONS (315mm)	
open width	580mm
open length	910mm » 920mm
open height	953mm
folded width	580mm
folded length	910mm » 978mm
folded height	607mm
frame wheelbase	370mm
REAR WHEEL	
Fixed	ALL MODELS
diameter	315mm,
CASTOR	
Solid	190mm

As a result of continual product improvement and development, RHealthcare reserve the right to change product specification without prior notice



R Healthcare, Building 2, Philips Campus,
Wellhall Road, Hamilton, ML3 9BZ
Tel: 0845 146 0600 Fax 01246 453 881
sales@rhealthcare.co.uk
www.rhealthcare.co.uk



ULTRA LIGHTWEIGHT ALUMINIUM CONSTRUCTION

WEIGHTS
12.3 KG

Dash Stowaway

*Easy folding, easy opening,
easy to carry, easy to store*



www.rhealthcare.co.uk

Dash Stowaway

*Easy folding, easy opening, easy to carry,
easy to store*

The Dash Stowaway with its quality lightweight aluminium frame represents great value to prescribers, carers and users alike. Backed by a 12 month parts guarantee.

Designed with help and support of users and distributors across the UK and around the world, the Dash Stowaway sees RHealthcare deliver a fantastic, lightweight, addition to the range of Dash wheelchairs, offering unrivalled value for money combined with high aesthetics and a unique, two way folding mechanism for easy lifting and storage.



Swing up arm rests
Aids easy side transfer



Flip up footplates
Length adjustable to fit the user comfortably



Simple quick flick catches
Folds fast and easy with no loose parts



Adjustable Brake Lever Length
Makes it easier to apply & release the brake



Ultra light and compact, the Dash Stowaway folding wheelchair allows the occupant to travel in style. Whether going on holiday, visiting the doctors, out for the day shopping or simply visiting friends, this highly adaptable, ultra compact, lightweight (only 12.3kgs) wheelchair is ideal for the occasional user both outdoors and indoors.

Easy folding, stores like a buggy unlike every other wheelchair meaning no more rear view restriction for car drivers, easily stores in the car boot or under the stairs. This chair is ready for use immediately as no parts need to be removed.

- **Stylish blue frame**
- **Built in tipping lever, ideal for kerbs**
- **315mm rear wheels with solid tyre**
- **Unique 2 way folding system**
- **17" Seat Width**
- **Nylon padded back upholstery for added comfort**
- **Lap belt option**
- **Fits 330 x 840 x 510 storage space, weighs less than 12.3kgs**
- **Low maintenance aluminium frame**
- **Maximum occupant weight limit 110kgs**

