



### Detachable Folding Powered Wheelchair

#### FreedomChair Model DE08L



With its clever detachable frame design, the DE08L separates into 2 sections with the heaviest part just 15.5kg. Yet it retains the classic sturdiness of FreedomChairs allowing active users up to 160kg to drive on uneven surfaces such as gravel, grass and cobblestones.

Suspension on front castors. Triple battery capacity up to 27 mile range.

#### Features

Seat Width: 18 inches (expandable)  
User Weight Limit: 160kg  
Drive Wheels: 12.5"

Motors: 2x 250W Continuous / 300W Peak  
Batteries: 1, 2 or 3 can be fitted with up to 9 miles/battery

Joystick controlled, FreedomChairs are versatile powered chairs for both outdoor use and, being highly manoeuvrable, excellent indoors too. They fold incredibly easily and quickly and fit into most car boots and take up very little space indoors. They are great travel chairs for days out, cruises etc and a battery test certificate is available for airline travel. With ISO, CE & FDA approval, the chairs are superbly made and finished. FreedomChairs can be fitted with multiple batteries, doubling or even tripling the range up to 27 miles!

---

#### Specifications

Weight: Seat section: 15.4kg. Rear wheel section: 13.8kg no batteries: 15.5kg with 1 battery  
Maximum weight capacity: 160kg  
Seat dimensions: 18" x 17.5"/45 x 43cm (expandable)  
Maximum speed: Up to 6 mph/9.5 kph  
Length: Open: 43.3"/110cm - Folded: 24"/61cm  
Width: Open & Folded: 23.9" / 61cm

Height: Open: 38.6"/98cm - Folded: 29.8"/75.6cm  
Gradient: 14 degrees  
Battery capacity: 24V 10AH per battery  
Drive wheels: 12.5" PU tyre  
Front wheels: 8" PU tyre with suspension  
Motors: 2x 250W Continuous / 300W Peak motors with electromagnetic braking