

TRAVELUX VENTURE Specification Sheet

Comfortable Captain seat with headrest

Width and height adjustable armrests

Adjustable seat height with elevating function

Fold-away one piece footplate



Pictured with Captain seat option*



Left or right hand joystick control option with seat height button



Easy fold-away one piece footplate



Fish-on seat option available

Available Frame Colour



Full Specifications

Total Length	102 cm	Turning Radius	63 cm
Total Width	64 cm	Maximum Speed	4 mph
Total Height	94 - 99 cm*	Range	26 km (16 miles)**
Seat Width	45 cm	Battery	2 x 12V / 35Ah
Seat Depth	40 cm	Charger	24V / 3Ah
Seat Height	53 - 69 cm*	Total Weight	71 kg
(From Ground)		(Including Battery)	
Backrest Height	48 cm	Maximum User Weight	135 kg (21 Stone)

KEY: * = Adjustable or optional ** = Varies depending on user weight, battery condition, weather and terrain conditions

www.vanosmedical.co.uk

