

# NEW FOR 2023!

## EXCEL QUEST POWERCHAIR

### Specification sheet

All new deluxe seat system!

Stadium style seat with a reclining backrest and super soft upholstery

90 – 115 degrees angle adjustable backrest

Quick release take apart chassis for easy transportation

Angle adjustable one piece footplate with increased plate size.

Available Colours 6 all new metallic colours:

Cherry Red,  
Sunburst Yellow,  
Solar Flare Orange,  
Envy Green,  
Paradise Blue,  
Satin Silver.



## Full Specification

Weight Capacity	136kg / 300lbs
Maximum Speed	4 mph
Range	16 miles
Batteries	2 x 18Ah
Length	82cm
Width	56cm
Height	105cm
Seat Height	54cm, 56.5cm, or 59 cm
Total Weight	57.6kg
Weight of heaviest part	17.5kg



Call us on  
+44 (0) 1977 681 400



[sales@vanosmedical.com](mailto:sales@vanosmedical.com)



[www.vanosmedical.co.uk](http://www.vanosmedical.co.uk)